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
Please read all information carefully.

Failure to properly follow the instructions may lead to serious surgical consequences.

Important: This package insert is designed to provide instructions for use of the ENSEAL[®] G2 Articulating Tissue Sealer. It is not a reference to surgical techniques.

ENSEAL[®] and I-BLADE[®] are trademarks of Ethicon Endo-Surgery.

**Instructions, Instructions, Gebrauchsanweisung, Istruzioni, Instruções, Instrucciones,
Gebruiksaanwijzing, Brugsvejledning, Ohje, Οδηγίες, Bruksanvisning, Instrukcja,
Utasítások, Návod k použití, Návod, 使用说明**

Indications

The ENSEAL® G2 Articulating Tissue Sealers are indicated for bipolar coagulation and mechanical transection of tissue during laparoscopic and open procedures.

The devices are bipolar electro-surgical instruments for use with the Generator G11 (GEN11). They are intended for use during open or laparoscopic, general and gynecological surgery to cut and seal vessels, and to cut, grasp and dissect tissue during surgery.

Indications for use include open and laparoscopic, general and gynecological surgical procedures (including urologic, thoracic, plastic and reconstructive, bowel resections, hysterectomies, cholecystectomies, gall bladder procedures, Nissen fundoplication, adhesiolysis, oophorectomies, etc.), or any procedure where vessel ligation (cutting and sealing), tissue grasping and dissection is performed. The devices can be used on vessels up to (and including) 7mm and tissue and/or vascular bundles as large as will fit in the jaws of the instruments.

The ENSEAL® G2 Articulating Tissue Sealers have not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use these devices for these procedures.

Contraindications

The efficacy of the ENSEAL® G2 Articulating Tissue Sealer for the indication of contraceptive tubal coagulation (permanent female sterilization) has not been evaluated and is unknown. The design of the ENSEAL® G2 Articulating Tissue Sealer is significantly different from bipolar designs that are marketed for the indication of contraceptive tubal coagulation. The design differences may affect the efficacy of the procedure and failure rates may not be comparable.

Device Description

The ENSEAL® G2 Articulating Tissue Sealer is a sterile, single patient use surgical instrument used to coagulate and transect vessels up to and including 7mm in diameter and tissue and/or vascular bundles. This device is for soft tissue only. The instrument consists of a grip housing assembly, a rotating and articulating shaft either 35cm or 45cm in length, a moveable jaw and an I-BLADE® in the jaws. The instrument shaft can be rotated to facilitate visualization and enable easy access to targeted tissue. The instrument shaft, normally in the straight position, can be articulated using the articulation wheel to gain additional access to tissue and facilitate additional angles across tissue and/or vascular bundles.

The jaws are in a normally opened position and can be partially or fully closed by squeezing the closing handle. The jaws have atraumatic teeth for grasping and holding targeted tissue when clamped. Bipolar energy is delivered to the clamped tissue when the energy activation button is pressed. Fully squeezing the closing handle advances the I-BLADE® the length of the jaws, and the targeted tissue is compressed, coagulated, and transected. The I-BLADE® cannot be advanced to transect tissue without first fully depressing the energy activation button. The power cord is permanently attached to the device and connects the instrument to the generator. The ENSEAL® G2 Articulating Tissue Sealer is designed for use exclusively with the Generator G11 (GEN11), packaged separately. Dimensions: Jaw Length, Straight: 20mm, Jaw Length, Curved: 19mm, Shaft Length: 35cm from Jaw Tip to Shaft Rotation Knob or 45cm from Jaw Tip to Shaft Rotation Knob, Shaft Diameter: 5mm.


Generator Audible Indicators

Refer to the following table for generator audible feedback tones.

Tone	Description	Action
Tone 1	Continuous steady beeping	Generator On: Device is active.
Tone 2	Continuous steady beeping at increased pitch	Intermediate: Tissue impedance threshold is reached.
Tone 3	One short beep	End of cycle: I-BLADE® fully advanced. Energy is stopped.

Illustration and Nomenclature (Illustration 1)

1. Jaws (curved or straight)
2. Atraumatic teeth
3. I-BLADE® knife (not visible)
4. Articulation joint
5. Shaft
6. Rotation knob
7. Articulation wheel
8. Energy activation button
9. Closing handle
10. Grip housing
11. Power cord

Transport and Storage Conditions

Temperature: -22° C to +60° C

Relative Humidity: 10% - 80%






How Supplied





The ENSEAL® G2 Articulating Tissue Sealer with integrated power cord is supplied sterile for single patient use. Discard after use.

Warnings and Precautions

- Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Examine the shipping carton and instrument for signs of shipping damage. Note any shortages, breakage, or apparent damage, retain the evidence, and notify Ethicon Endo-Surgery Customer Service immediately and replace with a new instrument. Do not use a damaged product.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised.
- Laparoscopic surgical procedures with the ENSEAL® G2 Articulating Tissue Sealer should be performed only by persons with adequate training and preparation. Consult the medical literature regarding techniques, complications, and hazards prior to any surgical procedure.
- Do not immerse the ENSEAL® G2 Articulating Tissue Sealer in liquid.
- Pools of conductive fluids (e.g. blood, saline) in direct contact with the instrument jaws may affect device performance. Remove pooled fluid from around the instrument jaws before activating the instrument.
- Verify compatibility with generators. Use device only with Ethicon Endo-Surgery Generator G11 (GEN11) software version X or later. Refer to the Generator G11 (GEN11) Operator's Manual to confirm software versions.
- The device is at high voltage whenever the generator is active. Activate this device only when the jaws are under direct visual observation to avoid unintended tissue coagulation or necrosis.
- Check if patient has a pacemaker or implanted cardioverter/defibrillator. Consult the manufacturers for information about the effects of RF energy on these devices.
- This instrument contains a surgically sharp edge. **THIS DEVICE CUTS!**
- The ENSEAL® G2 Articulating Tissue Sealer has not been shown to be effective for sealing large vessels without fully advancing the I-BLADE®. Do not attempt to seal large vessels by only clamping the tissue and applying energy.
- Do not use this device on vessels in excess of 7mm in diameter.

- Do not use excessive force on the closing handle to close the jaws; grasp only as much tissue as will fit between the jaws where the current will pass. Greater amounts of tissue require more closing handle force. Excessive force could damage the device.
- Avoid contact with metal and plastic when instrument is active.
- Use surgical gloves specifically designed by their manufacturer for electrosurgical procedures. The use of other gloves can result in burns or shocks.
- In case of system failure, ensure the availability of the appropriate back-up equipment relevant to the specific procedure. To articulate the device, turn the articulation wheel.
- During and following activation in tissue, the jaws and shaft may be hot. Avoid unintended contact with tissue, drapes, surgical gowns, or other unintended sites at all times.
- After using the instrument to seal and/or transect, always examine the tissue for hemostasis. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
- Due to the difficulty of visualizing obscured structures, proceed slowly and do not attempt to transect large masses of tissue in one activation. Avoid the division of large tissue and/or vascular bundles when using the instrument under these conditions.
- Products manufactured or distributed by companies not authorized by Ethicon Endo-Surgery may not be compatible with the ENSEAL® G2 Articulating Tissue Sealer instrument. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Discontinue use and replace device if excessive smoke or sparking is observed.
- Dispose of this device in a “Sharps” container. Dispose of all opened instruments whether used or unused.
- As with all energy sources (Electrosurgery, Laser, or Ultrasound), there are concerns about the carcinogenic and infectious potential of the by-products, such as tissue smoke plume and aerosols. Appropriate measures such as protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in laparoscopic or open procedures.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- To ensure effective sealing and transection of large vessels, always completely advance the I-BLADE® and wait for Tone 3.
- Do not stop squeezing the closing handle with the I-BLADE® only partially advanced.
- Transection without power application can lead to bleeding in vascular tissue.
- Do NOT immerse distal tip in liquid, including saline solution, to clean.
- Do not clean the instrument with abrasives.
- To avoid inadvertent user injury, do not activate the instrument while cleaning.
-  Prior to inserting or removing the device from the trocar, ensure that the device is in the straight, non-articulated position by confirming that the indicator arrow on the articulation wheel is aligned with the top of the instrument.
- Removing the device from a trocar while the device is articulated could result in unintentionally removing the trocar from the incision site, a loss of insufflation, or damage to the black shaft covering.
- If the black shaft covering is damaged, discontinue use.
- Prior to articulating the shaft, ensure that the jaw and the articulation joint are fully inserted through the trocar.
- The ENSEAL® G2 Articulating Tissue Sealer is not a passive articulation device. Do not force the articulation angle by pressing the shaft against a solid body. Use the articulation wheel to achieve the desired articulation angle.
- If the rotation force of the articulation wheel seems excessive, discontinue rotation of the wheel. Using excessive force could damage the device.
- Do not use the top jaw alone to perform tissue dissection.
- To reduce the possibility of an accidental activation, do not place your finger on the activation button until the jaws are in the desired position and activation is intended.
- This device is packaged and sterilized for single use only. Multiple patient use may compromise the device integrity or create a risk of contamination that, in turn, may result in patient injury or illness.

Instructions for Use	Illustrations & Notes
<p>1 Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field. Do not use a damaged product.</p> <p>2 Connect the instrument power cord to the generator.</p> <p>3 Before inserting the instrument through the trocar, make certain that the instrument is in the straight (non-articulated) position and the jaws are closed. (Illustration 2) Insert the instrument into the trocar before articulating the device.</p> <p>To rotate the device, turn the rotation knob in the desired direction. (Illustration 2a)</p> <p>To articulate the device, turn the articulation wheel. (Illustration 2b) The illustration provided shows the direction the shaft articulates when the rotation knob arrow is aligned on top.</p> <p>The articulation wheel may be turned even if the rotation knob indicator arrow is NOT aligned on top. The articulation wheel and rotation knob may be used independently.</p> <p>4 Position targeted tissue within open jaws. (Illustration 3).</p> <p>5 Squeeze the closing handle to grasp targeted tissue between the jaws (Illustration 4).</p> <p>6 Press the energy activation button to energize the jaws and begin coagulation of the targeted tissue. The generator gives audible feedback (Tone 1) as energy is delivered to the grasped tissue. Affirm tissue effects prior to advancing the I-BLADE® through tissue.</p>	<ul style="list-style-type: none"> NOTE: The instrument is straight when the indicator arrow on the articulation wheel is aligned with the top of the instrument.      <ul style="list-style-type: none"> NOTE: Tissue between the jaws may be further compressed by squeezing the closing handle to the first stop. Applying additional pressure at the first stop does not increase tissue compression.

Instructions for Use	Illustrations & Notes
<p>7 Complete the instrument firing cycle (Illustration 5).</p> <p>7a. Fully press the energy activation button to unlock the I-BLADE®.</p> <p>7b. While pressing the energy activation button, continue squeezing the closing handle to advance the I-BLADE® and increase tissue compression. NOTE: For larger vascular bundles, it is recommended to advance the I-BLADE® slowly (gradually).</p> <p>7c. The generator tone changes pitch (Tone 2) when the tissue impedance threshold is reached, and the I-BLADE® is ready to be fully advanced. It is recommended to maintain a consistent speed throughout I-BLADE® advancement.</p> <p>7d. Squeeze the closing handle until it stops against the grip housing. This fully advances the I-BLADE® (Illustration 5)</p> <p>The generator provides one short beep (Tone 3) to indicate the cycle is complete after the I-BLADE® is fully advanced and the tissue impedance threshold has been reached. Energy delivery automatically stops after Tone 3.</p>	<p>5</p>  <ul style="list-style-type: none"> NOTE: The instrument may be operated with either footswitch or handswitch control. NOTE: The I-BLADE® cannot be advanced to transect tissue without first fully pressing the energy activation button and then releasing the I-BLADE®. NOTE: The instrument can be wiped with a sterile moist gauze sponge to remove tissue, if necessary. If tissue is still visible in the jaws, use forceps to remove residue. (Illustration 6) <p>6</p>  <p>7</p> 
<p>8 After the cycle is complete, release the closing handle and the energy activation button. (Illustration 7) If the desired cut is longer than the stroke of the blade, transection can be accomplished by repeating Steps 4 - 7.</p>	<p>7</p> 
<p>9 Straighten the shaft and close the jaws before removing the device from the trocar. The instrument is straight when the indicator arrow on the articulation wheel is aligned with the top of the instrument.</p>	<p>8</p> 
<p>10 Disconnect the instrument from the generator.</p>	
<p>11 Dispose of the instrument in an appropriate Sharps container. Do not reuse instrument.</p>	
	<p>Caution: Do NOT immerse distal tip in liquid, including saline solution, to clean. (Illustration 8)</p>

This page intentionally left blank.

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ETHICON ENDO-SURGERY, LLC ©EES, LLC 2014
Ethicon Endo-Surgery (Europe) GmbH
Ethicon Endo-Surgery, Inc.